

Manufacturer: ChoicemedUSA, LLC.

Model: OxyWatch C316

Description: Fingertip Pulse Oximeter



Technical Specifications:

OLED Display

SPO2

Pulse Rate (PR)

SPO2 Measurement Range 70-99%

Resolution $\pm 1\%$

Accuracy 80-99% $\pm 2\%$, 70-79% $\pm 3\%$, < 69% Offers No Readings

Pulse Rate Range 30-235 BPM

Resolution ± 1 bpm

Accuracy ± 2 bpm

Battery Built-In Rechargeable Battery with Included AC/DC Adapter

Includes AC/DC Adapter, USB Cable and Software CD

One Year Manufacturer Limited Warranty

Reference: [Choicemed OxyWatch C316 Pulse Oximeter](#)