

Manufacturer: Contec Medical Systems

Model: CMS-50L

Description: Finger Pulse Oximeter



Technical Specifications:

Commonly Used in Hospitals, Other Medical Facilities, Homes, Fitness and Sports

Accurately Measures SPO2 and Pulse Rate (PR)

LED SPO2 and PR Display with Bar Graph

Low Battery Indicator

Automatic Shut Off When No Signal Detected

Pulse Rate Range 30-250BPM

Resolution ± 1 BPM

Accuracy $\pm 2\%$

SPO2 Range 35-99%

Resolution 1%

Accuracy $\pm 2\%$ 70-99%

Powered with 2 AAA Alkaline Batteries (30 Hours of Continuous Use)

Measures 66mm x 36mm x 33mm (Length x Width x Height)

Reference: Contec [CMS-50L Finger Pulse Oximeter](#)